

# REPROCESSOR CHEMICAL COMPATIBILITY


August 31, 2015

The ADVANTAGE PLUS<sup>®\*</sup>, ADVANTAGE PLUS<sup>®</sup> Pass-Thru<sup>\*\*</sup>, RAPIDAER<sup>®\*\*</sup>, ISA<sup>™</sup>, DSD EDGE<sup>®\*</sup>, DSD-201, SSD-102, CER OPTIMA<sup>®</sup> and CER series of endoscope reprocessors manufactured by MEDIVATORS have been tested and found to be compatible with the following high-level disinfectants (HLD) & detergents:

Product Name	HLD/Detergent	Reprocessor
RAPICIDE <sup>®</sup> PA Single Shot	HLD 5% peracetic acid	ADVANTAGE PLUS, ADVANTAGE PLUS Pass-Thru, DSD EDGE & RAPIDAER (UK Only) Reprocessors
RAPICIDE <sup>®</sup> Glut	HLD 2.5% acidic glutaraldehyde	DSD-201, SSD-102, CER OPTIMA & CER Series Reprocessors
RAPICIDE <sup>®</sup> OPA/28	HLD 0.575% <i>ortho</i> -Phthalaldehyde	DSD-201, SSD-102, CER OPTIMA & CER Series Reprocessors
RAPISPOR <sup>®</sup> PA	HLD 5% peracetic acid	RAPIDAER (Outside UK) Reprocessor
ADASPOR <sup>®</sup> /ADASPOR <sup>®</sup> Plus PA Single Shot	HLD 5% peracetic acid	ADVANTAGE PLUS & DSD EDGE Reprocessors
ADASPOR <sup>®</sup> /ADASPOR <sup>®</sup> Plus PA Ready-to-Use	HLD 5% peracetic acid	DSD-201 Reprocessor
RAPICIDE <sup>®</sup> PA Ready-to-Use	HLD 5% peracetic acid	DSD-201 & SSD-102 Reprocessors
ISASPOR <sup>®</sup> PA	HLD 5% peracetic acid	ISA Reprocessor
INTERCEPT <sup>®</sup> Detergent	Detergent	ADVANTAGE PLUS Reprocessor
INTERCEPT <sup>®</sup> Plus Detergent	Detergent	RAPIDAER, ADVANTAGE PLUS Pass-Thru Reprocessors

MEDIVATORS reserves the right to determine compatibility with our equipment. MEDIVATORS will not be responsible for damage to any equipment, including but not limited to MEDIVATORS reprocessors or flexible endoscopes of any brand, if this equipment has been exposed to chemicals not approved, endorsed or deemed compatible by MEDIVATORS.

- \* Note: ADVANTAGE PLUS and DSD EDGE Reprocessors have only been validated with INTERCEPT Detergent in the cleaning cycle.
- \*\* Note: RAPIDAER and ADVANTAGE PLUS Pass-Thru Reprocessors have only been validated with INTERCEPT PLUS Detergent in the cleaning cycle.



Craig Smith

Senior Vice President, Regulatory Affairs and Quality Assurance  
Corporate Compliance Officer  
Medivators Inc.

MEDIVATORS<sup>®</sup>, DSD EDGE<sup>®</sup>, ADVANTAGE PLUS<sup>®</sup>, RAPISPOR<sup>®</sup>, INTERCEPT<sup>®</sup>, and RAPICIDE<sup>®</sup> are registered trademarks of Medivators Inc.  
ADASPOR<sup>®</sup>, ISASPOR<sup>®</sup> and ISA<sup>™</sup> are registered trademarks and trademarks of Cantel Medical (Italy) S.r.l.  
RAPIDAER<sup>®</sup> is a registered trademark of Cantel Medical (UK) Limited.



Medivators Inc.  
14605 28th Avenue North  
Minneapolis, MN  
55447-4822 USA  
Tel: +1.763.553.3300  
Fax: +1.763.553.3387

Medivators BV  
Sourethweg 11  
6422 PC Heerlen  
The Netherlands  
Tel: +31.45.5.471.471  
Fax: +31.45.5.429.595

Cantel Medical Asia/Pacific Pte. Ltd.  
1A International Business Park  
#05-01 Singapore 609933  
Tel: +65.6227.9698  
Fax: +65.6225.6848

Cantel Medical Devices (China) Co. Ltd.  
Unit 804-805, Innov Tower Block A,  
Hongmei Road, Xuhui 200233 Shanghai  
Tel: +86 21 60161380  
Fax: +86 21 61210913

[www.medivators.com](http://www.medivators.com)