

# MEDIVATORS™

## Endoscope Drying Cabinet

## Endoscope Storage Cabinet



High efficiency,  
reliability and economy

Total microbiological  
security

Guaranteed dry  
scopes within 3 hours\*


Up to 31 days  
validated storage

# EVOLUTIONARY PROCESS, REVOLUTIONARY SOLUTION

For over 20 years, Cantel Medical's endoscope decontamination solutions have delivered unparalleled levels of efficiency, reliability and economy to over 300 endoscopy units worldwide.


It's therefore no wonder our MEDIVATORS™ Endoscope Drying & Storage Cabinet range has become the cabinet of choice, renowned for its reliability and unrivalled performance. In fact, the MEDIVATORS Endoscope Drying Cabinet is the only cabinet in the world fully compliant with EN-16442.

On average only 1 in 3 clean scopes are used on patients before reprocessing – an Endoscope Drying Cabinet halves that ratio, therefore paying for itself in just one year.




\*within the Endoscope Drying Cabinet

## Endoscope Drying Cabinet - Performance Data



## Endoscope Drying Cabinet - Temperature and Humidity Control



## FEATURES OF OUR ENDOSCOPE DRYING CABINET (EDC) AND ENDOSCOPE STORAGE CABINET (ESC) RANGE

- Dual filtered air maintains positive cabinet pressure
- HEPA filtered and dehumidified air delivered to each scope channel
- Clear toughened glass doors allow easy scope recognition
- Secure mounting for scope control section and lightguide plug
- All variations offer full traceability
- Extended storage time (up to 31 days) to cover holiday periods
- Four storage cabinet models available: 5 full size endoscopes, 10 full size endoscopes, 5 TOE probes, 10 ENT endoscopes
- Drying cabinet available for 10 full size endoscopes



## KEY USER BENEFITS – ENDOSCOPE STORAGE CABINET

- Up to 31 days validated storage
- Full traceability of scopes and operators
- Advanced fully Independent Monitoring System (IMS)


## KEY USER BENEFITS – ENDOSCOPE DRYING CABINET

In addition, the Endoscope Drying Cabinet offers the following features:

- The only drying cabinet fully compliant with EN-16442
- Dries scopes within 3 hours
- Temperature and humidity monitoring within the cabinet



## Product Specifications Cabinet 10T



### Specification

Model	Cabinet 10T
Dimensions	1300mm (w) x 610mm (d) x 2040mm (h)
Scope Capacity	10 scopes
Power Requirements	230v 50/60Hz
Lighting	2x low energy 15 watt lamps
Noise Level	< 58 dBA
Weight	450 kg

### Consumables

Description	Part number
10 Scope Drying Cabinet	FG-LA-EDC10
5 Scope Storage Cabinet	FG-LA-ESC5T
10 Scope Storage Cabinet	FG-LA-ESC10T
ESC ENT Cabinet	FG-LA-ESCENT
ESC TOE Cabinet	FG-LA-ESCTOE
Button Basket	LA-SBB
Main Override Key	GS-70201
Drip Mat	750002
Printer Paper	IS795732
Colonoscope Hanger	LA-ESCHA
Operator TAG	LA-ASGNT3
Scope TAG	LA-ASIST1

MEDIVATORS™ is a trademark of Medivators Inc.



Cantel Medical (UK) Limited  
 Wolseley House, Staffordshire Technology Park  
 Beaconside, Stafford, Staffordshire, ST18 0GA  
 United Kingdom  
 Tel: +44 (0) 1785.782420  
 Fax: +44 (0) 1785.782421

Medivators Inc.  
 14605 28th Avenue North  
 Minneapolis, MN  
 55447-4822 USA  
 Tel: +1.763.553.3300  
 Fax: +1.763.553.3387

Medivators BV  
 Sourethweg 11  
 6422 PC Heerlen  
 The Netherlands  
 Tel: +31.45.5.471.471  
 Fax: +31.45.5.429.595

Cantel Medical Asia/Pacific Pte. Ltd.  
 1A International Business Park  
 #05-01 Singapore 609933  
 Tel: +65.6227.9698  
 Fax: +65.6225.6848

Medivators Inc. Beijing Representative Office  
 Room 1801, Floor 18th, Tower A,  
 Beijing Marriott Hotel, Office Building  
 No. 7, Jianguomen South Avenue,  
 Dongcheng District, Beijing 100005 China  
 Tel: +86 10 65204039