

<i>Distribution:</i>	
International Distributors	✓
International Management	✓
US Clinical Specialists	✓
US Customer Service	✓
US Field Sales	✓
US Technical and Field Service	✓

Date: **October 21st, 2009**

RE: Deletion of Olympus® Channel Plug, Changes to the Air/Water/Gas Plug in Hook-up Assemblies and Creation of New Hook-ups

This Product Bulletin has been produced to provide notification that the Olympus Channel Plug (Olympus part number MH-944), and Air/Water/Gas Plug (Olympus part number MH-947) for Olympus 130 series and below endoscopes, will no longer be included with the following DSD hook-ups: DSD-110-HU0068, DSD-110-HU0069, DSD-110-HU0114 and DSD-110-HU0161. these hook-ups are now specifically intended for use with Olympus 140 series and above endoscopes.



The Olympus Channel Plug is required for manual cleaning of Olympus 140 series and above endoscopes, but has been removed from Medivators hook-up assemblies because it is provided by Olympus with every Olympus 140 series and above endoscope purchase. If necessary the Channel Plug can be ordered directly from Olympus, part number MH-944 or from Medivators, part number CLM-110-HU108.

Medivators has released four new hook-ups for reprocessing Olympus endoscopes 130 series and below. Olympus 130 series and below endoscopes will include the Air/Water/Gas Plug, which is required for reprocessing these endoscopes models. The new hook-up part numbers are DSD-110-HU0172, DSD-110-HU0173, DSD-110-HU0174, DSD-110-HU0175. These new hook-ups must also be used in-conjunction with the Channel Plug.

Due to the creation of the four new hook-ups for Olympus 130 series and below endoscopes, the following hook-ups should only be used for Olympus 140 series and above endoscopes: DSD-110-HU0068, DSD-110-HU0069, DSD-110-HU0114 and DSD-110-HU0161.

Please review the Hook-up Application Guide, PN 50097-111A for a cross-reference chart of endoscopes and applicable hook-ups. DFUs are available on the internet in the Resource Library at www.minntech.com/medivators.

If you have any questions related to this bulletin, please contact Medivators Technical Services Department at 1-800-444-4729.