

Certificate of Analysis

RAPICIDE[®] PA

Part A

Product Code: **ML02-0147**

Lot Number: 00016918

Expiration Date: 2019/11

Components	Specifications	Analysis
H ₂ O ₂	20.0 – 24.0% (w/w)	<u>22.0 %</u>
PAA	4.0 – 6.0 % (w/w)	<u>5.7 %</u>
Acetic Acid	8.0 – 10.0% (w/w)	<u>8.7 %</u>
pH, (as is)	0.5 – 1.1 (as is)	<u>0.6</u>


This lot has been reviewed and meets specified requirements.



Quality Analyst

11 May 2018

Analysis Date



Quality Assurance

11 MAY 18

Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

RAPICIDE[®] is a registered trademark of Medivators Inc.