



Certificate of Analysis

INTERCEPT[®] PLUS

Detergent & Cleaner


Part No. ML02-0145

Lot Number: 00026317

Expiration Date: 2019/09

Property	Test Method	Specification	Analysis
Appearance	Visual	Clear Non-viscous Liquid	Pass
pH (release) at 25°C	20002-154	11.5 – 12.5	12.1
Phosphate Level	20602-105	80.83 – 98.79 g/L	91.27 g/L

This lot has been reviewed and meets specified requirements.



Quality Analyst

05 Sept 2017
Analysis Date



Quality Assurance

06 Sept 2017
Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

INTERCEPT[®] is a registered trademark of Medivators Inc.