

# Certificate of Analysis

## **RAPICIDE<sup>®</sup> PA**

### **Part A**

Product Code: **ML02-0143**

Lot Number: 00032517

Expiration Date: 2019/04

Components	Specifications	Analysis
H <sub>2</sub> O <sub>2</sub>	20.0 – 24.0% (w/w)	<u>22.7 %</u>
PAA	4.0 – 6.0 % (w/w)	<u>5.5 %</u>
Acetic Acid	8.0 – 10.0% (w/w)	<u>8.7 %</u>
pH, (as is)	0.5 – 1.1 (as is)	<u>0.5</u>

This lot has been reviewed and meets specified requirements.

\_\_\_\_\_  
Quality Analyst

\_\_\_\_\_  
20 October 2017

\_\_\_\_\_  
Analysis Date

\_\_\_\_\_  
Quality Assurance

\_\_\_\_\_  
23-Oct-17

\_\_\_\_\_  
Release Date

**Medivators Inc.**  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)

RAPICIDE<sup>®</sup> is a registered trademark of Medivators Inc.