



MEDIVATORS
A Cantel Medical Company

Document
Number:

50092-648

Rev.
Level

E

TITLE

FORM, RAPICIDE OPA-28 TS

| REV | DESCRIPTION | Effective Date | |
|-----|--|----------------|-------|
| | | ORIGINATOR | ECO |
| D | ADD PASS/FAIL | SPW | 32460 |
| E | Update form to include test solutions, change info entry to lists/editable fields, correct spelling error. | PBS | 32975 |

Part Number ML02-0146

Description Rapicide OPA-28 Test Strips

Lot number Vial:039017 Kit:001017

Tested by Lyndee L. Heidecker

Date 23Jan2017

Note: Test form to be scanned and sent to incoming inspectors.

| Vial 1 | | Test Strip Color | |
|---------------|-------------|------------------|-------------|
| Concentration | Replicate 1 | Replicate 2 | Replicate 3 |
| 0.50 % | Pass | Pass | Pass |
| 0.35 % | Pass | Pass | Pass |
| Vial 2 | | Test Strip Color | |
| 0.50 % | Pass | Pass | Pass |
| 0.35 % | Pass | Pass | Pass |

| Visual Inspection | | Testing Method: Unaided Eye |
|-----------------------------|-------------|-----------------------------|
| | Pass / Fail | Comment |
| Cleanliness | Pass | Clean |
| Seal of caps | Pass | Sealed |
| Number of test strips -0/+3 | Pass | 50 Strips per vial |
| Adherence of test strip pad | Pass | Adhered |

| Testing Solutions | |
|-------------------|-----------------------|
| Concentration | Lot# / Exp Date |
| 0.50 % | CL16314-2 / 29Jan2017 |
| 0.35 % | CL16314-3 / 29Jan2017 |

CONFIDENTIAL AND PROPRIETARY INFORMATION OF MEDIVATORS INC. COPYING AND DISSEMINATION PROHIBITED, UNLESS AUTHORIZED BY MEDIVATORS. RETURN REQUIRED UPON DEMAND.