

Certificate of Conformance

Medline[®] OPA/28

Test Strips

Part Number:	MDS88OPASTP
Lot Number:	001115
Vial Lot Number:	029115
Expiration Date:	07/17

Quality Control Report:

Medivators certifies that the test strips indicate a positive or pass reaction with a solution mixed to a concentration of 0.50% \pm 0.05% OPA.

Medivators certifies that the test strips indicate a negative or fail reaction with a solution mixed to a concentration of 0.35% \pm 0.05% OPA.

The statistical analysis and approval test procedure is on file and available for review.

This lot has been reviewed and meets specified requirements.



Quality Assurance

1/18/2016

Release Date

Medivators Inc.
14605 28th Avenue N
Minneapolis, MN 55447
Tel: 800.444.4729
Fax: 763.553.3387
www.medivators.com

Medline[®] is a registered trademark of Medline.