



Quality Control Report and Certificate of Analysis

Formula: SB-130

Lot #: 001F08

	<u>Specification</u>	<u>Analysis</u>	<u>Deviation</u>
Sodium	<u>4500</u> mEq/L	<u>4476.5</u> mEq/L	<u>-0.52%</u>
Calcium	<u>0</u> mEq/L	<u>0.0</u> mEq/L	<u>0.00%</u>
Potassium	<u>0</u> mEq/L	<u>0.0</u> mEq/L	<u>0.00%</u>
Magnesium	<u>33.75</u> mEq/L	<u>33.3</u> mEq/L	<u>-1.19%</u>
Chloride	<u>4533.75</u> mEq/L	<u>4494.0</u> mEq/L	<u>-0.88%</u>
Dextrose	<u>9000</u> mg/dL	<u>9000.0</u> mg/dL	<u>0.00%</u>
Conductivity 25 C	<u>168.5</u> mS/cm	<u>166.9</u> mS/cm	<u>-0.95%</u>
Pyrogen Test	<u>Non-pyrogenic</u>	<u>Non-pyrogenic</u>	
pH Conc.	<u>1.8 - 3</u> pH Units	<u>1.99</u> pH units	
pH 1:34	<u>3 - 4</u> pH Units	<u>3.45</u> pH units	

Approved for Release:

20600-379/D

Deirdre
Quality Control

06-03-08
Date