



MEDIVATORS
A Cantel Medical Company

Certificate of Conformance

Medline[®] OPA/28

Test Strips

Part Number:	MDS88OPASTP
Lot Number:	002096
Vial Lot Number:	036096
Expiration Date:	2018-05-31

Quality Control Report:

Medivators certifies that the test strips indicate a positive or pass reaction with a solution mixed to a concentration of 0.50% $\pm$ 0.05% OPA.

Medivators certifies that the test strips indicate a negative or fail reaction with a solution mixed to a concentration of 0.35% $\pm$ 0.05% OPA.

The statistical analysis and approval test procedure is on file and available for review.

This lot has been reviewed and meets specified requirements.

Quality Assurance

9/21/2016

Release Date

Medivators Inc.
14605 28th Avenue N
Minneapolis, MN 55447
Tel: 800.444.4729
Fax: 763.553.3387
www.medivators.com

Medline Industries, Inc. is a registered trademark of Medline.