

## Certificate of Conformance

# Medline<sup>®</sup> OPA/28

## Test Strips

Part Number:	<b>MDS88OPASTP</b>
Lot Number:	003096A
Vial Lot Number:	036096A
Expiration Date:	2018-05-31

### Quality Control Report:

Medivators certifies that the test strips indicate a positive or pass reaction with a solution mixed to a concentration of 0.50%  $\pm$  0.05% OPA.

Medivators certifies that the test strips indicate a negative or fail reaction with a solution mixed to a concentration of 0.35%  $\pm$  0.05% OPA.

The statistical analysis and approval test procedure is on file and available for review.

This lot has been reviewed and meets specified requirements.



Quality Assurance

11/9/2016

Release Date

### Medivators Inc.

14605 28<sup>th</sup> Avenue N  
Minneapolis, MN 55447  
Tel: 800.444.4729  
Fax: 763.553.3387  
[www.medivators.com](http://www.medivators.com)

Medline Industries, Inc. is a registered trademark of Medline.