



Quality Control Report and Certificate of Analysis

Formula: SB-1003

Lot #: 003D07

	<u>Specification</u>	<u>Analysis</u>	<u>Deviation</u>
Sodium	<u>2946.4</u> mEq/L	<u>2908.5</u> mEq/L	<u>-1.29%</u>
Calcium	<u>128.9</u> mEq/L	<u>127.7</u> mEq/L	<u>-0.93%</u>
Potassium	<u>73.7</u> mEq/L	<u>72.6</u> mEq/L	<u>-1.46%</u>
Magnesium	<u>36.8</u> mEq/L	<u>36.3</u> mEq/L	<u>-1.36%</u>
Chloride	<u>3185.8</u> mEq/L	<u>3175.1</u> mEq/L	<u>-0.34%</u>
Dextrose	<u>7366</u> mg/dL	<u>7300.0</u> mg/dL	<u>-0.90%</u>
Conductivity 25 C	<u>151.5</u> mS/cm	<u>150.8</u> mS/cm	<u>-0.46%</u>
Pyrogen Test	<u>Non-pyrogenic</u>	<u>Non-pyrogenic</u>	
pH Conc.	<u>2 - 3</u> pH Units	<u>2.25</u> pH units	
pH 1:34	<u>3 - 4</u> pH Units	<u>3.55</u> pH units	

Approved for Release:

Jennifer Oestreich 4/3/07
Quality Control Date