



Quality Control Report and Certificate of Analysis

Formula: SB-1089

Lot #: 003F08

	<u>Specification</u>	<u>Analysis</u>	<u>Deviation</u>
Sodium	<u>3056.9</u> mEq/L	<u>3053.8</u> mEq/L	<u>-0.10%</u>
Calcium	<u>92.08</u> mEq/L	<u>91.6</u> mEq/L	<u>-0.58%</u>
Potassium	<u>147.32</u> mEq/L	<u>148.4</u> mEq/L	<u>0.73%</u>
Magnesium	<u>36.83</u> mEq/L	<u>36.5</u> mEq/L	<u>-1.03%</u>
Chloride	<u>3333.12</u> mEq/L	<u>3292.3</u> mEq/L	<u>-1.22%</u>
Dextrose	<u>7366</u> mg/dL	<u>7400.0</u> mg/dL	<u>0.46%</u>
Conductivity 25 C	<u>157.2</u> mS/cm	<u>155.9</u> mS/cm	<u>-0.83%</u>
Pyrogen Test	<u>Non-pyrogenic</u>	<u>Non-pyrogenic</u>	
pH Conc.	<u>2 - 3</u> pH Units	<u>2.1</u> pH units	
pH 1:34	<u>3 - 4</u> pH Units	<u>3.52</u> pH units	

Approved for Release:

20600-379/D


Quality Control

06-03-08
Date