



Quality Control Report and Certificate of Analysis

Formula: SB-1078

Lot #: 004D07

	<u>Specification</u>	<u>Analysis</u>	<u>Deviation</u>
Sodium	<u>2983.2</u> mEq/L	<u>2945.2</u> mEq/L	<u>-1.27%</u>
Calcium	<u>92.1</u> mEq/L	<u>90.4</u> mEq/L	<u>-1.85%</u>
Potassium	<u>36.8</u> mEq/L	<u>36.7</u> mEq/L	<u>-0.14%</u>
Magnesium	<u>36.8</u> mEq/L	<u>36.3</u> mEq/L	<u>-1.49%</u>
Chloride	<u>3149</u> mEq/L	<u>3136.6</u> mEq/L	<u>-0.39%</u>
Dextrose	<u>7366</u> mg/dL	<u>7300.0</u> mg/dL	<u>-0.90%</u>
Conductivity 25 C	<u>151</u> mS/cm	<u>148.9</u> mS/cm	<u>-1.39%</u>
Pyrogen Test	<u>Non-pyrogenic</u>	<u>Non-pyrogenic</u>	
pH Conc.	<u>2 - 3</u> pH Units	<u>2.24</u> pH units	
pH 1:34	<u>3 - 4</u> pH Units	<u>3.46</u> pH units	

Approved for Release:

Jennifer Oestreich
Quality Control

4/3/07
Date