



## Quality Control Report and Certificate of Analysis

Formula: SB-1014

Lot #: 004G07

|                   | <u>Specification</u>  | <u>Analysis</u>      | <u>Deviation</u> |
|-------------------|-----------------------|----------------------|------------------|
| Sodium            | <u>2983.2</u> mEq/L   | <u>3003.0</u> mEq/L  | <u>0.66%</u>     |
| Calcium           | <u>92.1</u> mEq/L     | <u>90.3</u> mEq/L    | <u>-2.01%</u>    |
| Potassium         | <u>0</u> mEq/L        | <u>0.0</u> mEq/L     | <u>0.00%</u>     |
| Magnesium         | <u>36.8</u> mEq/L     | <u>37.0</u> mEq/L    | <u>0.41%</u>     |
| Chloride          | <u>3112.1</u> mEq/L   | <u>3163.4</u> mEq/L  | <u>1.65%</u>     |
| Dextrose          | <u>7366</u> mg/dL     | <u>7200.0</u> mg/dL  | <u>-2.25%</u>    |
| Conductivity 25 C | <u>147.5</u> mS/cm    | <u>149.5</u> mS/cm   | <u>1.36%</u>     |
| Pyrogen Test      | <u>Non-pyrogenic</u>  | <u>Non-pyrogenic</u> |                  |
| pH Conc.          | <u>2 - 3</u> pH Units | <u>2.26</u> pH units |                  |
| pH 1:34           | <u>3 - 4</u> pH Units | <u>3.55</u> pH units |                  |

Approved for Release:

20600-379/D

Zachary Orr  
Quality Control

07-03-07  
Date