

## *Certificate of Conformance*

# RAPICIDE<sup>®</sup> OPA-28 Test Strips

**Part Number:**

**ML02-0136**

**Lot Number:**

005102A

**Expiration Date:**

2014/04

### **Quality Control Report:**

Minntech certifies that the test strips indicate a positive or pass reaction with a solution mixed to a concentration of  $50 \pm 10$  ppm Peracetic acid.

Minntech certifies that the test strips indicate a negative or fail reaction with a solution mixed to a concentration of  $5 \pm 1$  ppm Peracetic acid.

The statistical analysis and approval test procedure is on file and available for review.

This lot has been reviewed and meets specified requirements.



Quality Assurance

January 25, 2013

Release Date

Fax: 763.553.3387

[www.medivators.com](http://www.medivators.com)