




Quality Control Report and Certificate of Analysis

Formula: SB-1020

Lot #: 005F08

	<u>Specification</u>	<u>Analysis</u>	<u>Deviation</u>
Sodium	<u>2946.4</u> mEq/L	<u>2922.5</u> mEq/L	<u>-0.81%</u>
Calcium	<u>92.1</u> mEq/L	<u>90.1</u> mEq/L	<u>-2.23%</u>
Potassium	<u>73.7</u> mEq/L	<u>73.5</u> mEq/L	<u>-0.27%</u>
Magnesium	<u>36.8</u> mEq/L	<u>36.0</u> mEq/L	<u>-2.04%</u>
Chloride	<u>3149</u> mEq/L	<u>3137.8</u> mEq/L	<u>-0.36%</u>
Dextrose	<u>7366</u> mg/dL	<u>7400.0</u> mg/dL	<u>0.46%</u>
Conductivity 25 C	<u>148</u> mS/cm	<u>148.4</u> mS/cm	<u>0.27%</u>
Pyrogen Test	<u>Non-pyrogenic</u>	<u>Non-pyrogenic</u>	
pH Conc.	<u>2 - 3</u> pH Units	<u>2.03</u> pH units	
pH 1:34	<u>3 - 4</u> pH Units	<u>3.51</u> pH units	

Approved for Release:



Quality Control

06-03-08

Date