




Quality Control Report and Certificate of Analysis

Formula: SB-1003

Lot #: 005G08

	<u>Specification</u>	<u>Analysis</u>	<u>Deviation</u>
Sodium	<u>2946.4</u> mEq/L	<u>2922.5</u> mEq/L	<u>-0.81%</u>
Calcium	<u>128.9</u> mEq/L	<u>128.5</u> mEq/L	<u>-0.35%</u>
Potassium	<u>73.7</u> mEq/L	<u>73.3</u> mEq/L	<u>-0.51%</u>
Magnesium	<u>36.8</u> mEq/L	<u>36.2</u> mEq/L	<u>-1.77%</u>
Chloride	<u>3185.8</u> mEq/L	<u>3152.9</u> mEq/L	<u>-1.03%</u>
Dextrose	<u>7366</u> mg/dL	<u>7400.0</u> mg/dL	<u>0.46%</u>
Conductivity 25 C	<u>151.5</u> mS/cm	<u>149.2</u> mS/cm	<u>-1.52%</u>
Pyrogen Test	<u>Non-pyrogenic</u>	<u>Non-pyrogenic</u>	
pH Conc.	<u>2 - 3</u> pH Units	<u>2.06</u> pH units	
pH 1:34	<u>3 - 4</u> pH Units	<u>3.57</u> pH units	

Approved for Release:



Quality Control

07-02-08

Date