



Certificate of Analysis

RAPICIDE™ PA

RTU Part A

Lot Number: 007719

Expiration Date: 2020/08

| Components | Specifications | Analysis |
|-------------------------------|-------------------|--------------|
| H ₂ O ₂ | 20.0 – 24.0% | <u>23.3%</u> |
| PAA | 5.0% ± 1.0% | <u>5.5%</u> |
| Acetic Acid | 8.0 – 10.0% | <u>8.9%</u> |
| pH, (as is) | 0.5 - 1.1 (as is) | <u>0.5</u> |

This lot has been reviewed and meets specified requirements.

Quality Analyst

28-FEB-19

Analysis Date

Quality Assurance

04-MAR-19

Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

RAPICIDE™ is a trademark of Medivators Inc.