

**Certificate of Analysis**

**RAPICIDE® PA**

**Part A**

Product Code: **ML02-115**

Lot Number: **009819**


Expiration Date: **2020/09**

| Components                    | Specifications     | Analysis |
|-------------------------------|--------------------|----------|
| H <sub>2</sub> O <sub>2</sub> | 20.0 – 24.0% (w/w) | 23.4 %   |
| PAA                           | 4.0 – 6.0 % (w/w)  | 5.3 %    |
| Acetic Acid                   | 8.0 – 10.0% (w/w)  | 9.0 %    |
| pH, (as is)                   | 0.5 – 1.1 (as is)  | 0.6      |

This lot has been reviewed and meets specified requirements.

  
\_\_\_\_\_  
Quality Analyst

**12-MAR-19**  
\_\_\_\_\_  
Analysis Date

  
\_\_\_\_\_  
Quality Assurance

**14-MAR-19**  
\_\_\_\_\_  
Release Date

**Medivators Inc.**  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)

RAPICIDE® is a registered trademark of Medivators Inc.