



Certificate of Conformance

RAPICIDE® OPA-28 Test Strips

Part Number: **ML02-0136**

Lot Number: 011102

Expiration Date: 2014/04

Quality Control Report:

Minntech certifies that the test strips indicate a positive or pass reaction with a solution mixed to a concentration of 50 ± 10 ppm Peracetic acid.

Minntech certifies that the test strips indicate a negative or fail reaction with a solution mixed to a concentration of 5 ± 1 ppm Peracetic acid.

The statistical analysis and approval test procedure is on file and available for review.

This lot has been reviewed and meets specified requirements.

A handwritten signature in cursive script, reading "Diane Riskewitz".

Quality Assurance

January 8, 2013

Release Date

Fax: 763.553.3387
www.medivators.com