




Quality Control Report and Certificate of Analysis

Formula: SB-111

Lot #: 011A08

	<u>Specification</u>	<u>Analysis</u>	<u>Deviation</u>
Sodium	<u>4500</u> mEq/L	<u>4492.3</u> mEq/L	<u>-0.17%</u>
Calcium	<u>135</u> mEq/L	<u>131.1</u> mEq/L	<u>-2.89%</u>
Potassium	<u>90</u> mEq/L	<u>89.2</u> mEq/L	<u>-0.83%</u>
Magnesium	<u>33.75</u> mEq/L	<u>33.8</u> mEq/L	<u>0.00%</u>
Chloride	<u>4758.75</u> mEq/L	<u>4781.6</u> mEq/L	<u>0.48%</u>
Dextrose	<u>9000</u> mg/dL	<u>8900.0</u> mg/dL	<u>-1.11%</u>
Conductivity 25 C	<u>172</u> mS/cm	<u>164.3</u> mS/cm	<u>-4.48%</u>
Pyrogen Test	<u>Non-pyrogenic</u>	<u>Non-pyrogenic</u>	
pH Conc.	<u>1.8 - 3</u> pH Units	<u>2.01</u> pH units	
pH 1:34	<u>3 - 4</u> pH Units	<u>3.51</u> pH units	

Approved for Release:



Quality Control

01-03-08

Date