

# *Certificate of Conformance*

## **MaxiCide<sup>®</sup> OPA 28 Test Strips**

**Part Number:** **112-6977**

**Kit Lot Number:** 015087

**Vial Lot Number:** 045087

**Expiration Date:** 2019-08-31

### **Quality Control Report:**

MEDIVATORS certifies that the test strips indicate a positive or pass reaction with a solution mixed to a concentration of 0.50%  $\pm$  0.05% OPA.

MEDIVATORS certifies that the test strips indicate a negative or fail reaction with a solution mixed to a concentration of 0.35%  $\pm$  0.05% OPA.

The statistical analysis and approval test procedure is on file and available for review.

This lot has been reviewed and meets specified requirements.



*Quality Assurance*

9/6/2017

*Release Date*

### **Medivators Inc.**

14605 28<sup>th</sup> Avenue N  
Minneapolis, MN 55447

Tel: 800.444.4729

Fax: 763.553.3387

[www.medivators.com](http://www.medivators.com)