

Certificate of Analysis

RAPICIDE[®] PA

Part A

Product Code: **ML02-0158**

Lot Number: **015819**

Expiration Date: **2020/10**

| Components | Specifications | Analysis |
|-------------------------------|--------------------|---------------|
| H ₂ O ₂ | 20.0 – 24.0% (w/w) | 23.4 % |
| PAA | 4.0 – 6.0 % (w/w) | 5.3 % |
| Acetic Acid | 8.0 – 10.0% (w/w) | 9.0 % |
| pH, (as is) | 0.5 – 1.1 (as is) | 0.7 |


This lot has been reviewed and meets specified requirements.



Quality Analyst



Analysis Date



Quality Assurance



Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

RAPICIDE[®] is a registered trademark of Medivators Inc.

Certificate of Analysis

RAPICIDE® PA

Part B

Product Code: **ML02-0158**

Lot Number: **015919**

Expiration Date: **2020/10**

| Components | Specification | Analysis |
|---------------|---------------------|-------------|
| Sodium | 743.0 - 909.0 mEq/L | 760.0 mEq/L |
| Benzotriazole | 0.9 – 1.1% w/w | 0.9 %w/w |
| pH | 11.4 to 12.4 as is | 11.5 |


This lot has been reviewed and meets specified requirements.



Quality Analyst



Analysis Date



Quality Assurance



Release Date

Medivators Inc.
14605 28th Avenue North
Plymouth, Minnesota 55447
Telephone: (763) 553-3300
Fax (763) 553-3387

RAPICIDE® is a registered trademark of Medivators Inc.