

**Certificate of Analysis**

**RAPICIDE<sup>®</sup> PA**

**Part A**


Product Code: **ML 02-0115**

Lot Number: 015917


Expiration Date: 2018/11

Components	Specifications	Analysis
H <sub>2</sub> O <sub>2</sub>	20.0 – 24.0% (w/w)	<u>23.2 %</u>
PAA	4.0 – 6.0 % (w/w)	<u>5.5 %</u>
Acetic Acid	8.0 – 10.0% (w/w)	<u>8.9 %</u>
pH, (as is)	0.5 – 1.1 (as is)	<u>0.5</u>

This lot has been reviewed and meets specified requirements.

  
\_\_\_\_\_  
Quality Analyst

26 May 2017  
\_\_\_\_\_  
Analysis Date

  
\_\_\_\_\_  
Quality Assurance *PA*

29 May 2017  
\_\_\_\_\_  
Release Date

**Medivators Inc.**  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)

RAPICIDE<sup>®</sup> is a registered trademark of Medivators Inc.