



Certificate of Analysis

RAPICIDE® PA

RTU Part B

Lot Number: 018817
Expiration Date: 2018/12
Quality Control Report: 170615

Components	Specifications	Analysis
Na+	35.5 – 43.4 mEq/L	38.3
Benzotriazole	0.05% ± 0.005%	0.046
pH, (as is)	10.6 – 11.6 (as is)	11.52

This lot has been reviewed and meets specified requirements.

Quality Analyst

20 June 2017
Analysis Date

Quality Assurance

22 June 2017
Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com