



Certificate of Analysis
RAPICIDE™ PA
RTU Part B

Lot Number: 018818
Expiration Date: 2019/11

Components	Specifications	Analysis
Na+	35.5 – 43.4 mEq/L	<u>39.9mEq/L</u>
Benzotriazole	0.05% ± 0.005%	<u>0.053%</u>
pH, (as is)	10.6 – 11.6 (as is)	<u>11.6</u>

This lot has been reviewed and meets specified requirements.

Quality Analyst

25 May' 18
Analysis Date

Quality Assurance

25 May 18
Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

RAPICIDE™ is a trademark of Medivators Inc.