

Certificate of Analysis


RAPICIDE[®] PA

Part A

Product Code: **ML02-115**
Lot Number: 020718
Expiration Date: 2019/12

Components	Specifications	Analysis
H ₂ O ₂	20.0 – 24.0% (w/w)	<u>22.9 %</u>
PAA	4.0 – 6.0 % (w/w)	<u>5.5 %</u>
Acetic Acid	8.0 – 10.0% (w/w)	<u>8.8 %</u>
pH, (as is)	0.5 – 1.1 (as is)	<u>0.7</u>

This lot has been reviewed and meets specified requirements.



Quality Analyst

20 June 2018
Analysis Date



Quality Assurance

20 JUN 18
Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

RAPICIDE[®] is a registered trademark of Medivators Inc.