

**MINNTECH®**


A Cantel Medical Company


*Certificate of Conformance****RAPICIDE TEST STRIP*****Part Number:** ML02-0120**Lot Number:** 021474**Expiration Date:** 2013/10

I certify that the Rapicide test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the Rapicide test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.



Quality Analyst

Analysis Date

Quality Assurance

Release Date**MINNTECH CORPORATION**

14605 28th Avenue North

Plymouth, Minnesota 55447

Telephone: (763) 553-3300

Fax (763) 553-3387