



MINNTECH®
A Cantel Medical Company

Certificate of Conformance

RAPICIDE TEST STRIP

Part Number: ML02-0120


Lot Number: 021826

Expiration Date: 2013/11

I certify that the Rapicide test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the Rapicide test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.



Quality Analyst

01/10/12

Analysis Date



Quality Assurance

01/10/12

Release Date

MINNTECH CORPORATION
14605 28th Avenue North
Plymouth, Minnesota 55447
Telephone: (763) 553-3300
Fax (763) 553-3387

RECEIVING REPORT

Date: 01/10/12

Receiving Report No. 237976

V ALBERT BROWNE LTD CHANCERY HOU (105612)

E 190 WATERSIDE ROAD

Purchase Order No. Y13904

N HAMILTON IND PARK

Date of Receipt 01-10-12

D LEICESTER, . LE5 1QZ

O UK

Received at RECDOCK

R 011441162768636

Received by GDS

Received from ALBERT BROWNE LTD CHANCERY H

* THIS PRODUCT IS ON QA HOLD *

* UNTIL A GREEN ACCEPTED LABEL *

* IS PHYSICALLY PLACED ON THE *

* PRODUCT AND/OR THE PRODUCT IS *

* RELEASED IN THE ERP SYSTEM *

Carrier uti

Contact:

Freight Bill

REQUISITIONER

REQ NO. DEPARTMENT

PROJ. NO.

BUYER

Robert Weigel

PACKSLIP NO.

ITEM

PART #/WO/ACCT

DESCRIPTION

UM

QTY RCVD

DEL#

REV

SCHD

DATE

QTY REQD

QTY APPL

BAL DUE

167470

1 ML02-0120|01

Disp: I RAP GLUT TST EA

620

2 E

01-11-12

240

240

0

STRPS 2 VIAL ONE

3 E

01-12-12

1250

380

870

SHIPPER/BOX CONTAINS (2)

VIALS OF INDICATOR TEST

STRIPS, EACH CONTAINING

60 STRIPS.

Matl Cert:

Vis Sample Plan C=0 , AQL = 1.0 , Vis Sample Size = 34

Dim Sample Plan N/A , AQL = N/A

Lot Nbr : 022000 - 620

Special Notes: 1) Route to INSPECT
2) Q4 Receiving Inspection required for this item
3) Material Certification required for this item

Default Location: 01-J05A

Division Code MV - Endoscope Reprocessing