



MINNTECH
A Cantel Medical Company

Certificate of Conformance

RAPICIDE TEST STRIP

Part Number:	ML02-0120
Lot Number:	023760
Expiration Date:	2015-01

I certify that the Rapicide test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the Rapicide test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.

Quality Analyst

2/12/2013

Analysis Date

Quality Assurance

2/12/2013

Release Date

MINNTECH CORPORATION
14805 28th Avenue North
Plymouth, Minnesota 55447
Telephone: (763) 553-3300
Fax (763) 553-3387