

*Certificate of Conformance***RAPICIDE® OPA-28 Test Strips****Part Number:****73-OPA28B****Vial Lot Number:**

024065 kit lot number: 018065

Expiration Date:

02/17

Quality Control Report:

Minntech certifies that the test strips indicate a positive or pass reaction with a solution mixed to a concentration of $0.50\% \pm 0.05\%$ OPA.

Minntech certifies that the test strips indicate a negative or fail reaction with a solution mixed to a concentration of $0.35\% \pm 0.05\%$ OPA.

The statistical analysis and approval test procedure is on file and available for review.

This lot has been reviewed and meets specified requirements.


Quality Assurance7/24/2015
Release Date**MEDIVATORS Inc.**
14605 28th Avenue N
Minneapolis, MN 55447
Tel: 800.444.4729
Fax: 763.553.3387
www.medivators.com