

**MINNTECH®**

A Cantel Medical Company

*Certificate of Conformance***RAPICIDE TEST STRIP****Part Number:** ML02-0120**Lot Number:** 024465**Expiration Date:** 2015-06

I certify that the Rapicide test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the Rapicide test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.

Quality Analyst

7/24/2013

Analysis Date

Quality Assurance

7/24/2013

Release Date

MINNTECH CORPORATION

14605 28th Avenue North
Plymouth, Minnesota 55447
Telephone: (763) 553-3300
Fax (763) 553-3387

	DOCUMENTATION FORM	DOCUMENT NUMBER
		50092-591 REVISION: B PAGE 2 OF 2

Part Number MW02-0120
 Description Appt cide Test strips
 Lot number 024465
 Tested by S. Walsh
 Date 7-24-13

VIAL/TEST LOCATION: <u>Sample</u>		TEST STRIP COLOR	
Concentration	Replicate 1 Pass/Fail	Replicate 2 Pass/Fail	Replicate 3 Pass/Fail
2.5 %	<u>Pass</u>	<u>Pass</u>	<u>Pass</u>
1.3 %	<u>Pass</u>	<u>Pass</u>	<u>Pass</u>

VIAL/TEST LOCATION:		TEST STRIP COLOR	
Concentration	Replicate 1 Pass/Fail	Replicate 2 Pass/Fail	Replicate 3 Pass/Fail
2.5 %			
1.3 %			<u>N/A</u>

TEST NO 2 – 5 Visual Inspection

TESTING METHOD : Unaided Eye

	Pass / Fail	Comment
Cleanliness	<u>Pass</u>	<u>clean</u>
Seal of caps	<u>Pass</u>	<u>Sealed</u>
Number of test strips	<u>Pass</u>	<u>60 Test strips</u>
Adherence of test strip path	<u>Pass</u>	<u>Adhered</u>