



MEDIVATORS
A Cantel Medical Company

Certificate of Conformance

RAPICIDE® Test Strip

Part Number: ML02-0120

Lot Number: 025669

Expiration Date: 2016-04

I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.

Quality Analyst

06/18/2014

Analysis Date

Quality Assurance

06/18/2014

Release Date

Medivators Inc.

14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com