



Certificate of Conformance

RAPICIDE® Test Strip

Part Number: ML02-0120

Lot Number: 026322

Expiration Date: 2016-10

I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.
I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.

Quality Analyst

20Nov2014

Analysis Date

Quality Assurance

20Nov2014

Release Date

Medivators Inc.

14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

**MEDIVATORS**

A Cantel Medical Company

Document
Number:**50092-591**Rev.
Level**D****TITLE****FORM, RAPICIDE GLUT TS**

REV	DESCRIPTION	Effective Date	ORIGINATOR	ECO
C	REVISE TO MATCH INSPECTION PROCEDURE		DR	31243
D	ADD PASS/FAIL		SPW	32460

Part Number ML02-0120**Description** Rapicide Glutaraldehyde Test Strips**Lot number** 026322 exp 2016-10**Tested by** MTV**Date** 19 Nov 14

Note: Test form to be scanned and sent to incoming inspectors.

vial 1	TEST STRIP COLOR		
Concentration	Replicate 1	Replicate 2	Replicate 3
2.5 %	Pass / Fail	Pass / Fail	Pass / Fail
1.3 %	Pass / Fail	Pass / Fail	Pass / Fail

vial 2	TEST STRIP COLOR		
Concentration	Replicate 1	Replicate 2	Replicate 3
2.5 %	Pass / Fail	Pass / Fail	Pass / Fail
1.3 %	Pass / Fail	Pass / Fail	Pass / Fail

TESTING METHOD: Unaided Eye

	Pass / Fail	Comment
Cleanliness	Pass / Fail	Clean
Seal of caps	Pass / Fail	Sealed
Number of test strips -0/+3	Pass / Fail	60 strips per vial
Adherence of test strip path	Pass / Fail	Adhered