



Certificate of Conformance

RAPICIDE® Test Strip

Part Number: ML02-0120


Lot Number: 026585

Expiration Date: 2016-11

I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.



Quality Analyst

12Jan2015

Analysis Date



Quality Assurance

12Jan2015

Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

**MEDIVATORS**

A Cantel Medical Company

Document
Number:**50092-591**Rev.
Level**E****TITLE****FORM, RAPICIDE GLUT TS**

		Effective Date	
REV	DESCRIPTION	ORIGINATOR	ECO
D	ADD PASS/FAIL	SPW	32460
E	Update form to include test solutions, change info entry to lists/editable fields, correct spelling error, include kit and vial lot fields.	PBS	32975

Part Number ML02-0120**Description** RapicideTest Strips**Lot number** Vial:026585 Kit:026585**Tested by** Nick Votrobeck**Date** 12Jan2015

Note: Test form to be scanned and sent to incoming inspectors.

Vial 1

Concentration	Replicate 1	Replicate 2	Replicate 3
2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass

Vial 2

2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass

Visual Inspection**Testing Method: Unaided Eye**

	Pass / Fail	Comment
Cleanliness	Pass	Clean
Seal of caps	Pass	Sealed
Number of test strips -0/+3	Pass	60/vial
Adherence of test strip pad	Pass	Adhered

Testing Solutions

Concentration	Lot# / Exp
2.5 %	760479 / Exp 23Jan2015
1.3 %	CL14357-3 / Exp 23Jan2015