



Certificate of Analysis

RAPICIDE[®] PA

RTU Part A

Lot Number: 028617

Expiration Date: 2019/03

Quality Control Report: 170893

Components	Specifications	Analysis
H ₂ O ₂	20.0 – 24.0%	<u>23.5</u>
PAA	5.0% ± 1.0%	<u>5.1</u>
Acetic Acid	8.0 – 10.0%	<u>9.1</u>
pH, (as is)	0.5 - 1.1 (as is)	<u>1.00</u>

This lot has been reviewed and meets specified requirements.

Quality Analyst

22 Sept 2017
Analysis Date

Quality Assurance

25 SEP 17
Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com