

# Certificate of Analysis

## **RAPICIDE® PA**

### **Part A**



Product Code: **ML02-115**  
Lot Number: **029019**  
Expiration Date: **2021/01**

Components	Specifications	Analysis
H <sub>2</sub> O <sub>2</sub>	20.0 – 24.0% (w/w)	<u>23.3 %</u>
PAA	4.0 – 6.0 % (w/w)	<u>5.1 %</u>
Acetic Acid	8.0 – 10.0% (w/w)	<u>9.1 %</u>
pH, (as is)	0.5 – 1.1 (as is)	<u>0.6</u>

This lot has been reviewed and meets specified requirements.

  
\_\_\_\_\_  
Quality Analyst

**24-LUG-19**  
\_\_\_\_\_  
Analysis Date

 **QP**   
\_\_\_\_\_  
Quality Assurance

**24/07/19**  
\_\_\_\_\_  
Release Date

**Medivators Inc.**  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)

RAPICIDE® is a registered trademark of Medivators Inc.