

**Certificate of Analysis**

**RAPICIDE® PA**

**Part A**

Product Code: ML02-115

Lot Number: 029219

Expiration Date: 2021/01

Components	Specifications	Analysis
H <sub>2</sub> O <sub>2</sub>	20.0 – 24.0% (w/w)	<u>22.6 %</u>
PAA	4.0 – 6.0 % (w/w)	<u>5.2 %</u>
Acetic Acid	8.0 – 10.0% (w/w)	<u>9.1 %</u>
pH, (as is)	0.5 – 1.1 (as is)	<u>0.6</u>

This lot has been reviewed and meets specified requirements.

FDC  
Quality Analyst

24-LUG-19  
Analysis Date

AP  
Quality Assurance

24/02/19  
Release Date

**Medivators Inc.**  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)

RAPICIDE® is a registered trademark of Medivators Inc.