



DOCUMENTATION  
STANDARD OPERATING  
PROCEDURE

DOCUMENT  
NUMBER

20100-716

REVISION: D  
PAGE 7 OF 7

VI. Quality Control report and certificate

Formula: BC-1-R  
Lot #: 02x07

Concentrate Specification

Sodium	<u>1187.7</u> mEq/l $\pm$ 2%
Chloride	<u>402.6</u> mEq/l $\pm$ 2%
Bicarbonate	<u>785.1</u> mEq/l $\pm$ 2%
pH	<u>7.0 – 8.3</u> pH units
Conductivity (@ 25°C)	<u>67</u> mS/cm $\pm$ 5%

When mixed using the correct equipment and dilution with AAMI certified water Formula BC-1-R,  
Lot # 02x07 meets the above specifications.

Approved for release:

Quality Control

10-11-07

Date