



# Certificate of Conformance

## **RAPICIDE® Test Strip**

**Part Number:** ML02-0120

**Lot Number:** 030290

**Expiration Date:** 2019-03-29

I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.

Quality Analyst

01MAY2017

Analysis Date

Quality Assurance

01MAY2017

Release Date

**Medivators Inc.**

14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)





**MEDIVATORS**  
A Cantel Medical Company

Document  
Number:

**50092-591**

Rev.  
Level

**E**

**TITLE**

**FORM, RAPICIDE GLUT TS**

**Effective Date**

**REV**

**DESCRIPTION**

**ORIGINATOR**

**ECO**

D

ADD PASS/FAIL

SPW

32460

E

Update form to include test solutions, change info entry to lists/editable fields, correct spelling error, include kit and vial lot fields.

PBS

32975

**Part Number** ML02-0120

**Description** RapicideTest Strips

**Lot number** Vial:030290 Kit:030290

**Tested by** Erik A. Grotte

**Date** 01MAY2017

Note: Test form to be scanned and sent to incoming inspectors.

**Vial 1**

Concentration	Replicate 1	Replicate 2	Replicate 3
2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass

**Vial 2**

2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass

**Visual Inspection**

**Testing Method: Unaided Eye**

	Pass / Fail	Comment
Cleanliness	Pass	CLEAN
Seal of caps	Pass	SEALED
Number of test strips -0/+3	Pass	60 STRIPS PER VIAL
Adherence of test strip pad	Pass	ADHERED

**Testing Solutions**

Concentration	Lot# / Exp
2.5 %	852591 / 2018-09-28
1.3 %	CL17121-1 / 01JUN2017