



A Cantel Medical Company

Certificate of Conformance

Part Number:

ML02-0120

Lot Number:

030894

Expiration Date:

2019-08-24

I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.


Quality Analyst

28Sep2017

Analysis Date

Huber *Cont*
Quality Assurance

28Sep2017

Release Date

Medivators Inc.

14605 28th Avenue North

Minneapolis, MN 55447 USA

Telephone: (763) 553-3300

Fax (763) 553-3387

www.medivators.com