



**MEDIVATORS**  
A Cantel Medical Company

## Certificate of Conformance

# **RAPICIDE<sup>®</sup> Test Strip**

**Part Number:** ML02-0120

**Lot Number:** 031406

**Expiration Date:** 2019-09-13

I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.

Quality Analyst

29Nov2017

Analysis Date

Quality Assurance

29Nov2017

Release Date

**Medivators Inc.**

14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)