



**Certificate of Analysis**

**RAPICIDE™ PA**

**RTU Part A**

Lot Number: 031619

Expiration Date: 2021/02

| Components                    | Specifications    | Analysis     |
|-------------------------------|-------------------|--------------|
| H <sub>2</sub> O <sub>2</sub> | 20.0 – 24.0%      | <u>23.5%</u> |
| PAA                           | 5.0% ± 1.0%       | <u>5.3%</u>  |
| Acetic Acid                   | 8.0 – 10.0%       | <u>9.0%</u>  |
| pH, (as is)                   | 0.5 - 1.1 (as is) | <u>0.6</u>   |

This lot has been reviewed and meets specified requirements.

  
\_\_\_\_\_  
Quality Analyst

21-A60-19  
\_\_\_\_\_  
Analysis Date

  
\_\_\_\_\_  
Quality Assurance

21/08/19  
\_\_\_\_\_  
Release Date

**Medivators Inc.**  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)

RAPICIDE™ is a trademark of Medivators Inc.