



Certificate of Conformance

RAPICIDE® Test Strip

Part Number: ML02-0120


Lot Number: 033133

Expiration Date: 2020-10-30

I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.



Quality Analyst

28Dec2018

Analysis Date



Quality Assurance

28Dec2018

Release Date

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

**MEDIVATORS**

A Cantel Medical Company

Document
Number:**50092-591**Rev.
Level**E****TITLE****FORM, RAPICIDE GLUT TS**

REV	DESCRIPTION	Effective Date	ORIGINATOR	ECO
D	ADD PASS/FAIL		SPW	32460
E	Update form to include test solutions, change info entry to lists/editable fields, correct spelling error, include kit and vial lot fields.		PBS	32975

Part Number ML02-0120**Description** RapicideTest Strips**Lot number** Vial:033133 Kit:033133**Tested by** Tong Vang**Date** 28Dec2018

Note: Test form to be scanned and sent to incoming inspectors.

Vial 1

Concentration	Replicate 1	Replicate 2	Replicate 3
2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass

Vial 2

2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass

Visual Inspection**Testing Method: Unaided Eye**

	Pass / Fail	Comment
Cleanliness	Pass	Clean
Seal of caps	Pass	Sealed
Number of test strips -0/+3	Pass	Vial1: 62 Strips ; Vial2: 60 Strips
Adherence of test strip pad	Pass	Adhered

Testing Solutions

Concentration	Lot# / Exp
2.5 %	942059 / 05Jan2019
1.3 %	CL18339-2 / 05Jan2019