



# Certificate of Conformance

## **RAPICIDE® Test Strip**

**Part Number:** ML02-0120


**Lot Number:** 033134

**Expiration Date:** 2020-10-30


I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.

  
\_\_\_\_\_  
Quality Analyst

28Dec2018  
\_\_\_\_\_  
Analysis Date

  
\_\_\_\_\_  
Quality Assurance

28Dec2018  
\_\_\_\_\_  
Release Date

**Medivators Inc.**  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)



**MEDIVATORS**  
A Cantel Medical Company

Document  
Number:

**50092-591**

Rev.  
Level

**E**

**TITLE**

**FORM, RAPICIDE GLUT TS**

REV	DESCRIPTION	Effective Date	ORIGINATOR	ECO
D	ADD PASS/FAIL		SPW	32460
E	Update form to include test solutions, change info entry to lists/editable fields, correct spelling error, include kit and vial lot fields.		PBS	32975

**Part Number** ML02-0120

**Description** RapicideTest Strips

**Lot number** Vial:033134 Kit:033134

**Tested by** Tong Vang

**Date** 28Dec2018

Note: Test form to be scanned and sent to incoming inspectors.

Vial 1			
Concentration	Replicate 1	Replicate 2	Replicate 3
2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass
Vial 2			
2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass

Visual Inspection		Testing Method: Unaided Eye
	Pass / Fail	Comment
Cleanliness	Pass	Clean
Seal of caps	Pass	Sealed
Number of test strips -0/+3	Pass	60 Strips per Vial
Adherence of test strip pad	Pass	Adhered

Testing Solutions	
Concentration	Lot# / Exp
2.5 %	942059 / 05Jan2019
1.3 %	CL18339-2 / 05Jan2019