



MEDIVATORS
A Cantel Medical Company

Certificate of Conformance
RAPICIDE® Test Strip

Part Number: ML02-0120

Lot Number: 033135

Expiration Date: 2020-10-31

I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.

Quality Analyst

14Feb2019

Analysis Date

Quality Assurance

14Feb2019

Release Date

Medivators Inc.

14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com



MEDIVATORS
A Cantel Medical Company

Document
Number:

50092-591

Rev.
Level

E

TITLE

FORM, RAPICIDE GLUT TS

REV	DESCRIPTION	Effective Date	ORIGINATOR	ECO
D	ADD PASS/FAIL		SPW	32460
E	Update form to include test solutions, change info entry to lists/editable fields, correct spelling error, include kit and vial lot fields.		PBS	32975

Part Number ML02-0120

Description RapicideTest Strips

Lot number Vial:033135 Kit:033135

Tested by Nick Votrobeck

Date 14Feb2019

Note: Test form to be scanned and sent to incoming inspectors.

Vial 1			
Concentration	Replicate 1	Replicate 2	Replicate 3
2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass
Vial 2			
2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass

Visual Inspection		Testing Method: Unaided Eye
	Pass / Fail	Comment
Cleanliness	Pass	Clean
Seal of caps	Pass	Sealed
Number of test strips -0/+3	Pass	Vial 1: 62 strips; Vial 2: 60 strips
Adherence of test strip pad	Pass	Adhered

Testing Solutions	
Concentration	Lot# / Exp
2.5 %	400080 / 14Mar2019
1.3 %	CL19045-2 / 14Mar2019

[Signature] 14Feb2019

CONFIDENTIAL AND PROPRIETARY INFORMATION OF MEDIVATORS INC. COPYING AND DISSEMINATION PROHIBITED, UNLESS AUTHORIZED BY MEDIVATORS. RETURN REQUIRED UPON DEMAND.