



Certificate of Conformance

RAPICIDE[®] Test Strip

Part Number: ML02-0120

Lot Number: 033136

Expiration Date: 2020-10-31

I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.

Quality Analyst

14Feb2019

Analysis Date

Quality Assurance

14Feb2019

Release Date

Medivators Inc.

14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com



MEDIVATORS
A Cantel Medical Company

Document
Number:

50092-591

Rev.
Level

E

TITLE

FORM, RAPICIDE GLUT TS

| REV | DESCRIPTION | Effective Date | ECO |
|-----|---|----------------|-------|
| D | ADD PASS/FAIL | SPW | 32460 |
| E | Update form to include test solutions, change info entry to lists/editable fields, correct spelling error, include kit and vial lot fields. | PBS | 32975 |

Part Number ML02-0120

Description RapicideTest Strips

Lot number Vial:033136 Kit:033136

Tested by Nick Votrobeck

Date 14Feb2019

Note: Test form to be scanned and sent to incoming inspectors.

| Vial 1 | | | |
|---------------|-------------|-------------|-------------|
| Concentration | Replicate 1 | Replicate 2 | Replicate 3 |
| 2.5 % | Pass | Pass | Pass |
| 1.3 % | Pass | Pass | Pass |
| Vial 2 | | | |
| 2.5 % | Pass | Pass | Pass |
| 1.3 % | Pass | Pass | Pass |

| Visual Inspection | | Testing Method: Unaided Eye |
|-----------------------------|-------------|-----------------------------|
| | Pass / Fail | Comment |
| Cleanliness | Pass | Clean |
| Seal of caps | Pass | Sealed |
| Number of test strips -0/+3 | Pass | 60 strips per vial |
| Adherence of test strip pad | Pass | Adhered |

| Testing Solutions | |
|-------------------|-----------------------|
| Concentration | Lot# / Exp |
| 2.5 % | 400080 / 14Mar2019 |
| 1.3 % | CL19045-2 / 14Mar2019 |

Nick Votrobeck 14 Feb 2019