



**Certificate of Analysis**

**RAPICIDE<sup>®</sup> PA**

**RTU Part A**

Lot Number: 033217

Expiration Date: 2019/04

Quality Control Report: 171006

Components	Specifications	Analysis
H <sub>2</sub> O <sub>2</sub>	20.0 – 24.0%	<u>22.3%</u>
PAA	5.0% ± 1.0%	<u>5.8%</u>
Acetic Acid	8.0 – 10.0%	<u>8.6%</u>
pH, (as is)	0.5 - 1.1 (as is)	<u>1.0</u>

This lot has been reviewed and meets specified requirements.

\_\_\_\_\_  
Quality Analyst

\_\_\_\_\_  
19 October 2017  
Analysis Date

\_\_\_\_\_  
Quality Assurance

\_\_\_\_\_  
23-OCT-17  
Release Date

**Medivators Inc.**  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)