



Certificate of Conformance

RAPICIDE® Test Strip

Part Number: ML02-0120

Lot Number: 034104

Expiration Date: 2021-04-02

I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.

The statistical analysis and approved test procedure is on file and available for review.

Quality Analyst

14May2019

Analysis Date

Quality Assurance

14May2019

Release Date

Medivators Inc.

14605 28th Avenue North
Minneapolis, MN 55447 USA
Telephone: (763) 553-3300
Fax (763) 553-3387
www.medivators.com

**MEDIVATORS**

A Cantel Medical Company

Document
Number:**50092-591**Rev.
Level**E****TITLE****FORM, RAPICIDE GLUT TS**

REV	DESCRIPTION	Effective Date	ECO
D	ADD PASS/FAIL	SPW	32460
E	Update form to include test solutions, change info entry to lists/editable fields, correct spelling error, include kit and vial lot fields.	PBS	32975

Part Number ML02-0120**Description** RapicideTest Strips**Lot number** Vial:034104 Kit:034104**Tested by** Nick Votrobeck**Date** 14May2019

Note: Test form to be scanned and sent to incoming inspectors.

Vial 1			
Concentration	Replicate 1	Replicate 2	Replicate 3
2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass
Vial 2			
2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass

Visual Inspection		Testing Method: Unaided Eye
	Pass / Fail	Comment
Cleanliness	Pass	Clean
Seal of caps	Pass	Sealed
Number of test strips -0/+3	Pass	Vial 1: 62, Vial 2: 60
Adherence of test strip pad	Pass	Adhered

Testing Solutions	
Concentration	Lot# / Exp
2.5 %	415730 / 14Jun2019
1.3 %	CL19134-1 / 14Jun2019

14 May 2019

CONFIDENTIAL AND PROPRIETARY INFORMATION OF MEDIVATORS INC. COPYING AND DISSEMINATION PROHIBITED, UNLESS AUTHORIZED BY MEDIVATORS. RETURN REQUIRED UPON DEMAND.